


Mtub02/RD105 Analysis----Beijing Family Identification

Principle

TB strains belonging to the Beijing family possess a specific deletion called RD105 by Tsolaki *et al.* (J. Clin. Microbiol. 2005. Vol. 43, p3185-3191) and that we independently observed using primers tagging the Mtub02 VNTR.



By combining two PCRs with primers either flanking the deletion Mtub02Beij (Mub02F + Mtub02BeR) or primers inside the deletion Mtub02 (Mtub02F + Mtub02R) it is possible to quickly separate isolates into Beijing-type and non-Beijing type.

On the figure above the samples 1, 7, 8, and 9 and the control strain H37Rv are non-Beijing type strains. The samples 2, 3, 4, 5, 6, 10, 11, 12, 13, and 14 are Beijing-type strains.

1. PCR amplification

1) Sample DNA

Use 2µl of DNA extract (5ng/µl) per amplification reaction.

2) Primers

(1) Mtub02

L): 5'-CGTGACAGTTGGGTGTTTA-3'

R): 5'-TTCGTTCAAGAACTCCAAGG-3'

(2) Mtub02Be

Mtub02RevBeijing CGTCGTTTTCTGCCGATACT

3) Preparation of working master mix.

(1) Preparation for the working master mix for 100 reactions

Reagent	Volume(μ l)	*100Rns(μ l)
10XBuffer(Including MgCl ₂)	1.5	150
dNTP	1.5	150
5MBetain	3.0	300
Primers(Left and Right mixture)	1.5	100
Taq E (5U/ μ l)	0.05	5
Distilled Water	5.5	600
Total	13.0	1300

(2) Dispense 13 μ l of working master mix to each PCR reaction tube.

(3) Add 2 μ l of DNA samples(5ng/ μ l) to each reaction vial as given below.

4) Thermal cycling protocol

The cycling protocol is as follows:

94°C 5min

94°C 30sec

62°C 30sec 35cycles

72°C 45sec

72°C 10min

20°C until the gel is loaded for electrophoresis

2. Detection and Analysis of PCR Products by Agarose Gel Electrophoresis(AGE).